

Med-Aire 8" Alternating Pressure and Low Air Loss Mattress System



INSET C

- Active mattress replacement system provides both alternating pressure and low air loss to optimize pressure redistribution, shear/friction reduction, and microclimate control
- Designed to prevent, treat and heal pressure ulcers in the home or long term care setting
- Control unit offers 10 minute cycles when alternating and can also be used in static mode
- Adjustable patient weight settings allow for optimal immersion and patient comfort & compliance
- Audio and Visual alarm alerts when air pressure is low
- System includes 20 individual air cells to maximize immersion and envelopment with pillow feature for head stability & comfort
- 8" deep air cells are easily removed and replaced with air leakage plugs
- Removable Nylon stretch cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- CPR valve allows for rapid deflation
- Tubing quickly and easily disconnects with EZ Lock/Quick Release design

PRODUCTS	Description
14027	Alternating Pressure/Low Air Loss Mattress System, E0277, 1/ea

SPECIFICATIONS	Description
Control Unit Airflow LPM	8 LPM
Control Unit Cycle Time	10 Min.
Control Unit Type	Compressor
Control Unit Alarms	Visual / Audible
Mattress Cover Material	Fluid-Resistant Nylon
Weight Capacity	350 lbs.
Warranty	1 Year Limited

ACCESSORIES	Description
14027M	Mattress Only, 1/ea
14027P	Pump Only, 1/ea
14027-CELL	Replacement Bladder for 14027 and 14028N, 1/ea
14027-CELL 1	Ventilated Bladder, 1/ea
14427	Electric Pump (Sold Separately) (Replaces Mason # 8151-EP), 1/ea
14027c	Cover for 14027, 1/ea

PARTS	Description
14027-BASE	Replacement Base for 14027, 1/ea
14027c	Cover for 14027, 1/ea
14027-CELL	Replacement Bladder for 14027 and 14028N, 1/ea
14027-cell1	Repl Bladder w/Lsr Holes, 1/ea

PARTS	Description
14028v	One Way Valve For Low Air Loss Matt, 1/ea
14012	CPR Plug ForLow AirLoss Matt, 1/ea
14028t	Tubing, 1/ea
14001T	Beige Tubing for 14001 and 14001E, 1/ea
14027-fftc	Connector for 14027, 1/ea
14001l	"L" Connector, 1/ea
14027-fmtc	Bladder Quick Release Connector, 1/ea
14027-fbs	Stopper for 14027, 1/ea
14028h	Hose, 1/ea
14027P	Pump Only, 1/ea
14028s	Pump socket, 1/ea
14007or-3	"O" Ring, 1/ea
14007	EZ Lock Q/R Connector, 1/ea
140pdk	Adjustment knob, 1/ea
14027-fk	Adjustment Knob, 1/ea
14027-fph	Hanging hook, 1/ea
14029f	Intake filter, 1/ea
14530pf	Intake filter, 1/ea
14530fc	Filter Cover, 1/ea
14530f	Fuse (1 amp 250V), 1/ea
14001-l	"L" Connector, 1/ea