

FX3 FREE BALANCE SELF STANDING FOLDING CANE



Item #	Product Description	HCPCS	UOM
GVFX310GY-6	FX3 Free Balance Folding Cane, Matte Grey	E0100	6/cs
GVFX310BK-6	FX3 Free Balance Folding Cane, Black	E0100	6/cs

FEATURES

- Cane base with 3 feet, pivots in any direction providing added security and traction
- 3 point base allow the cane to stand on its own
- Tip flexes and swivels to maintain contact with the ground when the cane is at an angle, providing superior balance
- Conveniently folds for easy storage in seconds
- Soft grip handle relieves stress and pressure on fingers and wrist
- Height adjusts for use by most individuals
- Tip is removable without tools, for easy replacement
- Weight Capacity 300 lbs

Specifications

Handle Height	32.5" x 39.5"
Number of Height Adjustments	8
Weight Capacity	300 lbs.
Warranty	Limited Lifetime